



At all times, even if you are listed below as being eligible to engage in a BioCharger session, ensure that you listen to your body. At no time should you feel any discomfort or pain, and if you do, we suggest you cease the session.

### Not Eligible to Engage the BioCharger.

- Pregnant Woman.
- Anyone with an electrical brain implant.
- Anyone with an artificial implant in the head of any kind.

### Eligible to Engage with Preventative Measures.

Steps to ensure safety and positive experience.  
Children younger than 12 years:

Age (years)	Minimum distance from BioCharger
8-12	1 meter (± 3 feet)
6-8	1.25 meter (± 4 feet)
4-6	1.5 meter (± 5 feet)
1-4	2 meter (± 6 feet)

#### Photo-Sensitive Individuals:

Epileptics, photo-convulsive individuals and anyone that is photo sensitive in any way should take one of the following preventative measures:

- Use an opaque eye mask to block all light to the eyes.
- Turn their back to the BioCharger, preventing light from reaching their eyes.
- Listen to your body and stop use if feeling uncomfortable.

### Eligible to Engage the BioCharger at a distance.

- Engage BioCharger at minimum distance of 2 meters (± 6 feet)
- People with pace makers or other implants for heart.
- Non-removable insulin pumps (removable pumps must be removed prior to starting session).
- Any other electrical implant.
- Any electrical device if not a removable device we recommend a 6 foot distance and consulting your physician prior to use..

### Protect the BioCharger

- Only use supplied stylus to operate screen in order to prevent electrical damage to screen.
- Ensure stylus is clean and dry.
- Don't touch BioCharger or screen while running a recipe.
- Don't place any liquids on or around BioCharger.
- Ensure vents on rear are clear of debris.
- Ensure BioCharger housing is clean and free of dust and grime.
- Don't operate if liquid has been spilled.
- Follow installation and operation guides at all times.

If using the BioCharger for the first time it is recommended to engage the session at a distance of 2 meters (±6 feet) and then move closer as you feel more comfortable and don't experience any discomfort.

It is important to note that the BioCharger is not a medical device and should not be considered a substitute for consulting with your health care professional. If you have any questions, please consult with your health care professional of choice.



#biochargerNG



#biochargerNG



#BioChargerNG



<http://biocharger.com>